

HIALUROM TENDON

INSTRUCTIONS FOR USE



HIALUROM Tendon

Sodium hyaluronate 40 mg/2 ml, sterile solution for injection in pre-filled syringe

For peritendinous or intrasheath injection.

HIALUROM Tendon is a viscoelastic, isotonic, 40 mg/2 ml solution of sodium hyaluronate in a buffered physiological solution comprised of sodium phosphate, sodium chloride, and mannitol at a pH of 7 - 8. The sodium hyaluronate is purified from the bacterial fermentation of a *Streptococcus* strain.

Sterile by moist heat.

Indications:

For the treatment of pain and restricted mobility in tendon disorders.

Contraindications:

- HIALUROM Tendon is contraindicated:
- for patients with a history of hypersensitivity (allergy) to one of its constituents
- for patients with a history of hypersensitivity to gram positive bacterial products
- for patients with an active infection or skin disease near or around the injection site
- for patients who have received another injectable medicine near or around the injection site

Potential adverse events:

Local secondary phenomena such as pain, feeling of heat, bruising, redness and swelling may occur following treatment with HIALUROM Tendon.

Dosage and administration:

HIALUROM Tendon should be administered only by medical specialists trained in peritendinous or intrasheath injection technique.

Intrasheath injection:

In tendons with a sheath, inject HIALUROM Tendon into the tendon sheath in the affected area. Peritendinous injection:

In tendons without a sheath, inject the solution along the affected tendon, but not in the tendon.

Inject HIALUROM Tendon around the affected tendon or into the affected tendon sheath once a week for a total of 2 injections. Several tendons may be treated at the same time. Repeat treatments may be administered as required.

The content of the HIALUROM Tendon pre-filled syringe is sterile.

Take the pre-filled syringe out of the pack, unscrew the cap, attach a suitable needle (e.g., 25- to 27-gauge) and secure by turning slightly.

The final needle selection for any procedure is determined by the physician. There are several factors, which need to be considered in choosing the size (gauge) and length of the needle to be used for peritendinous or intrasheath injection, including the anatomy of the region, distance between the skin and tendon, and patient characteristics (weight, age). Ultrasound guidance is recommended during injection.

Remove any air bubble, if present, before injection.

Interactions:

No information on the incompatibility of HIALUROM Tendon with other medications administered to tendons is available to date.

Precautions:

Caution should be exercised in patients with known hypersensitivity to medicinal products. The general precautions for peritendinous and intrasheath injections should be observed. HIALUROM Tendon should be instilled accurately into the tendon sheath or around the affected tendon. The ultrasound guidance is recommended during injection.

Avoid nerve lesions and injections into blood vessels!

The procedure should be avoided in patients with known systemic bleeding disorders, or in patients with history of vasovagal reaction or syncope.

As no clinical evidence is available on the use of sodium hyaluronate in children and adolescents, pregnant and lactating women, treatment with HIALUROM Tendon is not recommended in these cases.

Patients who have experienced any complications in the days after injection should contact the physician immediately.

Warnings:

Single use device! Each pre-filled syringe of HIALUROM Tendon 40 mg/2 ml, solution for injection is intended to be used once only for a single patient.

Do not use if the pre-filled syringe or the blister are damaged. Any solution not used immediately after opening must be discarded. Otherwise, the sterility is no longer guaranteed.

The used needles and syringes must be discarded after injection and should not be kept for other administrations. Do not re-use! Reuse of syringes or needles already used can lead to the transmission of infectious agents (including HIV and hepatitis).

Do not resuscitate, as this may damage or alter the product.

Strict aseptic administration technique must be followed to minimize infections at the injection site.

Injection site must be properly disinfected (70% alcohol or with another disinfectant). Disinfectants containing quaternary ammonium salts should not be used for skin preparation as hyaluronic acid can precipitate under such conditions.

Do not use if the product is yellow, discoloured or contains precipitates. Keep out of the reach of children!

HIALUROM Tendon is a medical device. To be used by a physician only.

Characteristics and mode of action:

A tendon is a strong structure of fibrous connective tissue designed to transmit forces from muscle to bone and resist load during muscle contraction. Tendons may be surrounded by different structures: fibrous bands, synovial sheaths, peritendon sheaths, tendon bursa. Overuse or inappropriate biomechanical stress may cause inflammation and/or degenerative changes of the tendon, leading to pain and loss of function. Lubricating the tendon could reduce pain, improve tendon function and reduce the potential for adhesions.

Because of its lubricating and viscoelastic properties, HIALUROM Tendon promotes tendon gliding and the physiological repair process. In addition, due to its macromolecular meshwork HIALUROM Tendon reduces the free passage of inflammatory cells and molecules.

HIALUROM Tendon is a transparent solution of natural and highly purified sodium hyaluronate obtained by fermentation and is devoid of animal protein. HIALUROM Tendon also contains mannitol, a free radical scavenger, which helps to stabilise the chains of sodium hyaluronate.

Presentation:

HIALUROM Tendon - One pre-filled syringe of 40 mg/2 ml.

Storage

Do not freeze. The product is stored at temperatures below 25°C, in original package. Do not use after the expiry date indicated on the box!

Last revision date: June 2016

Manufacturer:

S.C. ROMPHARM COMPANY S.R.L.
1A Eroilor Street, Olopeni, Ilfov, 075100, Romania

Distributor:

LIBYTEC PHARMACEUTICAL S.A.
24, Voulagmenis Avenue, 167 77 Ellinika, Greece.
PHONE: 210 96 09 960, FAX 210 96 38 438



Explanation of Symbols

	Consult Instructions for Use
	Batch code
	Use by date
	Do not resuscitate
	Do not re-use
	Sterilised using steam or dry heat
	Upper limit of temperature
	Do not use if package is damaged
	Manufacturer
	Product conform with requirements in the European Medical Devices Directive
	Number of Notified Body



HIALUROM TENDON



HIALUROM Tendon

Sodium hyaluronate 40 mg/2 ml, sterile solution for injection in pre-filled syringe.

For peritendinous or intrasheath injection.

To HIALUROM Tendon is a viscoelastic, isotonic, 40 mg/2 ml solution of sodium hyaluronate in a buffered physiological solution comprised of sodium phosphate, sodium chloride, and mannitol at a pH of 7 - 8. The sodium hyaluronate is purified from the bacterial fermentation of a *Streptococcus* strain.

Aseptic packaging by moist heat.

Indications:

For the treatment of pain and restricted mobility in tendon disorders.

Contraindications:

- HIALUROM Tendon is contraindicated:
- for patients with a history of hypersensitivity (allergy) to one of its constituents
- for patients with a history of hypersensitivity to gram positive bacterial products
- for patients with an active infection or skin disease near or around the injection site
- for patients who have received another injectable medicine near or around the injection site

Potential adverse events:

Local secondary phenomena such as pain, feeling of heat, bruising, redness and swelling may occur following treatment with HIALUROM Tendon.

Dosage and administration:

HIALUROM Tendon should be administered only by medical specialists trained in peritendinous or intrasheath injection technique.

Intrasheath injection:

In tendons with a sheath, inject HIALUROM Tendon into the tendon sheath in the affected area. Peritendinous injection:

In tendons without a sheath, inject the solution along the affected tendon, but not in the tendon.

Inject HIALUROM Tendon around the affected tendon or into the affected tendon sheath once a week for a total of 2 injections. Several tendons may be treated at the same time. Repeat treatments may be administered as required.

The content of the HIALUROM Tendon pre-filled syringe is sterile.

Take the pre-filled syringe out of the pack, unscrew the cap, attach a suitable needle (e.g., 25- to 27-gauge) and secure by turning slightly.

The final needle selection for any procedure is determined by the physician. There are several factors, which need to be considered in choosing the size (gauge) and length of the needle to be used for peritendinous or intrasheath injection, including the anatomy of the region, distance between the skin and tendon, and patient characteristics (weight, age). Ultrasound guidance is recommended during injection.

Remove any air bubble, if present, before injection.

Interactions:

No information on the incompatibility of HIALUROM Tendon with other medications administered to tendons is available to date.

Precautions:

Caution should be exercised in patients with known hypersensitivity to medicinal products. The general precautions for peritendinous and intrasheath injections should be observed. HIALUROM Tendon should be instilled accurately into the tendon sheath or around the affected tendon. The ultrasound guidance is recommended during injection.

Avoid nerve lesions and injections into blood vessels!

The procedure should be avoided in patients with known systemic bleeding disorders, or in patients with history of vasovagal reaction or syncope.

As no clinical evidence is available on the use of sodium hyaluronate in children and adolescents, pregnant and lactating women, treatment with HIALUROM Tendon is not recommended in these cases.

Warnings:

Single use device! Each pre-filled syringe of HIALUROM Tendon 40 mg/2 ml, solution for injection is intended to be used once only for a single patient.

Do not use if the pre-filled syringe or the blister are damaged. Any solution not used immediately after opening must be discarded. Otherwise, the sterility is no longer guaranteed.

The used needles and syringes must be discarded after injection and should not be kept for other administrations. Do not re-use! Reuse of syringes or needles already used can lead to the transmission of infectious agents (including HIV and hepatitis).

Do not resuscitate, as this may damage or alter the product.

Strict aseptic administration technique must be followed to minimize infections at the injection site.

Injection site must be properly disinfected (70% alcohol or with another disinfectant). Disinfectants containing quaternary ammonium salts should not be used for skin preparation as hyaluronic acid can precipitate under such conditions.

Do not use if the product is yellow, discoloured or contains precipitates. Keep out of the reach of children!

HIALUROM Tendon is a medical device. To be used by a physician only.

Characteristics and mode of action:

A tendon is a strong structure of fibrous connective tissue designed to transmit forces from muscle to bone and resist load during muscle contraction. Tendons may be surrounded by different structures: fibrous bands, synovial sheaths, peritendon sheaths, tendon bursa. Overuse or inappropriate biomechanical stress may cause inflammation and/or degenerative changes of the tendon, leading to pain and loss of function. Lubricating the tendon could reduce pain, improve tendon function and reduce the potential for adhesions.

Because of its lubricating and viscoelastic properties, HIALUROM Tendon promotes tendon gliding and the physiological repair process. In addition, due to its macromolecular meshwork HIALUROM Tendon reduces the free passage of inflammatory cells and molecules.

HIALUROM Tendon is a transparent solution of natural and highly purified sodium hyaluronate obtained by fermentation and is devoid of animal protein. HIALUROM Tendon also contains mannitol, a free radical scavenger, which helps to stabilise the chains of sodium hyaluronate.

Presentation:

HIALUROM Tendon - One pre-filled syringe of 40 mg/2 ml.

Storage

Do not freeze. The product is stored at temperatures below 25°C, in original package. Do not use after the expiry date indicated on the box!

Manufacturer:

S.C. ROMPHARM COMPANY S.R.L.
1A Eroilor Street, Olopeni, Ilfov, 075100, Romania

Distributor:

LIBYTEC PHARMACEUTICAL S.A.
A. Boudouris Avenue, 167 77 Ellinika, Elladika
TEL: 210 96 09 960, FAX 210 96 38 438



Explanation of Symbols

	Consult Instructions for Use
	Batch code
	Use by date
	Do not resuscitate
	Do not re-use
	Sterilised using steam or dry heat
	Upper limit of temperature
	Do not use if package is damaged
	Manufacturer
	Product conform with requirements in the European Medical Devices Directive
	Number of Notified Body



fiber direction

Elaborat	Verificat și aprobat
Departament achiziții	Întreținere
Sef RA/Manager de contract	

Font: Times New Roman 9 pt.; Spatiu text 3,516 mm

TECHNICAL DETAILS	CLIENT	ROMPHARM	DETINATOR APP	ROMPHARM	
	PRODUCT	Hialurom Tendon		COD INTERN	3023257D10
	DIMENSION	140 x 600 mm (Folding in the middle)		TARA	—
	COLOR	BLACK		PH-CODE	2524
	MATERIAL	op medical 40g/m ²			
	BRAILLE	—			
	CLISEU	—			
VARNISH	—				

ABC - IMPEX	DIECUT		CYLINDER	
	DTP		PROOF READER	
	DATA	21.06.2018		
APPROVAL	APPROVED			
	NOT APPROVED, RESEND			
	Pentru o vizualizare corecta functia "overprint preview" trebuie sa fie activata! For accurate visualisation please activate the overprint preview.			
	SIGNATURE			